

MTS3000 - RINGCORE

Automatic System for Residual Stress

Measurement by Ring-Core



MTS3000 - RingCore

Automatic System for Residual Stress

Measurement by Ring-Core

Main features

- Automatic control of the feed motion of the drilling tool by stepper motor
- Automatic control of the drilling
- Digital camera & laser beams for positioning at the centre of the strain gage
- Hollow shaft with a strain gage rosette special quick installation tool
- Measurement of residual stress variation with depth (from 0.25mm to 5.0mm)
- Choice of residual stress calculation methods
- · Automatic zero-depth detection
- High measuring accuracy
- Dedicated software for control of the measurement (RSM-RC) and back-calculation of residual stress (EVAL-RC)





Residual stresses and measurement by the ring core method


The ring-core method is a partially-destructive test method for measuring the residual stresses present in a mechanical component and consists in machining an annular groove around a special strain gage rosette.

It can be used for determining uniform and non-uniform residual stresses in depth, mainly on large-size forgings or castings: it consists in machining a core (14mm internal diameter / 18mm external diameter / 5mm depth) around a special strain gage rosette with a three overlapped grids.


Compared with all other methods for residual stress measurements, the ring-core method has some distinctive features: it enables the profile of residual stresses in a material up to a depth of about 5 mm. One of the main advantages of the ring-core method is also the high sensitivity (3 to 4 times higher than hole drilling), which reduces the influence of the error in the acquisition process. The influence of core-rosette eccentricity error is also lower than hole-drilling as the eccentricity errors compensate each other.

The MTS3000 – RingCore automatic system facilitates the measurement procedure and allows to obtain more accurate results than with manual systems.

A couple of laser beams and an high resolution camera allows the centering of the strain gage with the drilling axis. Moreover, for a quick installation of the strain gage, a dedicated installation tool can be used.












MTS3000 – RINGCORE

Automatic System for Residual Stress Measurement by Ring-Core

The MTS3000-RingCore system for measuring residual stresses

The MTS3000-RingCore system is composed by:

- Mechanical device (drilling unit for hollow milling cutter, digital camera & laser beams for centering)
- Electronic devices (control of the advancing and drilling, zero setting surface)
- Dedicated software (software for the automatic control and the back-calculation of residual stress)
- Digital strain gage amplifier (automatic acquisition of the strain gage rosette data)



MTS3000 - RingCore: drilling control & acquisition software (RSM-RC)

The MTS3000-RingCore system is automatically controlled by the RSM-RC software.

The software allows the control of the drilling unit, the setting of the measurement (strain gage selection, amplifier) and allows also the automatic acquisition of the strain values from the rosette during the measurement.







The software allows also the acquisition of the digital camera installed on the mechanical body of the device: this camera is used for the centering the strain gage with the drilling axis thanks to the presence of two laser beam that are able to project two perpendicular lines directly on the tested workpiece. These lines can be aligned with the strain gage reference positioning lines.


MTS3000-RingCore: residual stress evaluation software (EVAL-RC)

Acquired data can be processed using the SINT EVAL-RC analysis software.

The EVAL-RC software allows the calculation of residual stress using the Differential method, the Incremental method and the Integral method. Moreover, the EVAL-RC software gives the possibility to make graphs of residual stress in the depth, calculating the stress in a custom direction or plotting a Mohr circle of residual stress.









Contact us

Via delle Calandre, 63 - 50041 Calenzano (FI) ITALY

Tel: +39.055.8826302

Fax: +39.055.8826303

www.sintechnology.com

info@sintechnology.com

VAT n. 04185870484

Recognitions

SINT Technology's test laboratory is accredited to standard ISO/IEC 17025:2005 by the Italian accreditation body **ACCREDIA with certificate no. 0910-L**





Certification of conformity to the requirements of standard

UNI EN ISO 9001

